

MYS-502N4

MECHANICAL HOSPITAL BED
WITH 2 CRANKS



The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font. A registered trademark symbol (®) is located at the top right of the word.

www.meyosis.com


MYS-502N4



GENERAL SPECIFICATIONS

- Vascular, legrest and knee movement positions are controlled by 2 cranks manually.
- Removable Polypropylene (PP) plastic head and foot ends with colour panels.
- Split PP plastic side rails with soft drop release mechanism and angle indicator.
- 4 independent and detachable ABS mattress bases.
- Antibacterial electrostatic powder coated steel frame, open chasis.
- 125 mm, 360 degree swivel wheels with cross brake
- Accessory adaptors on the corners, 1 height adjustable IV pole.
- Urine bag holders and belt hooks on both sides.
- Protective plastic bumpers on corners
- Alternative colour options for side rails
- 26 Density foam mattress with antibacterial, waterproof and zippered cover.

DIMENSIONS



TECHNICAL SPECIFICATIONS



OPERATION ANGLES

