

MYS-503N

MANUAL HOSPITAL BED WITH
TRENDLENBURG



Meyosis[®]

www.meyosis.com

MYS-503N



GENERAL SPECIFICATIONS

- Backrest, legrest and knee movements, vascular positions are controlled by 2 crank
- Manuel trendelenburg
- Removable Polypropylene (PP) plastic head and foot ends with colour panels
- Foldable aluminum side rails
- 4 independent and detachable ABS lying sections
- Antibacterial electrostatic powder coated steel frame, open chasis,
 - Ø125 mm, 360° swivel wheels with cross brake system
- Accesory adaptors on the corners, 1 height adjustable IV Pole
- Urine bag holders and belt hooks on both sides
- Protective plastic bumpers on corners
- 32 Density foam mattress with antibacterial, waterproof and zippered cover

OPERATION ANGLES



Vascular
0° - 30°



Knee Movement



Back Movement
0° - 70°




Leg Movement
0° - 30°



Trendelenburg
0° - 17°

DIMENSIONS



TECHNICAL SPECIFICATIONS



100 kg



Max. 200 kg



Manual