

MYS-511N

BABY CRIB



Meyosis[®]

www.meyosis.com

MYS-511N

GENERAL SPECIFICATIONS



— Chrome-plated metal frame

— 360-degree rotatable swivel castors with cross brakes

— High-quality plexiglass lying surface

— Waterproof, zippered mattress

— Trendelenburg position

— Wire basket underneath

— Weight: 20kg
Maximum capacity: 40kg

DIMENSIONS

- Height Min: 22 cm
 Max: 26 cm
- Basket Insides: 39,5 cm
 Outsides: 45 cm
- Length Outsides: 76 cm
 Insides: 69cm
- Profile 25*2 mm, 20*30