

MYS-520N

PEDIATRIC HOSPITAL BED WITH
2 MOTORS



The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font. A registered trademark symbol (®) is located at the top right of the word.

www.meyosis.com


MYS-520N

GENERAL SPECIFICATIONS

- Bedside, calf movement is adjusted with 2 dc motors.
- Positions are controlled by hand control.
- Footrest movement is 10 steps manually with destomat gear.
- Electrostatic powder coated metal body.
- Full length gradual working metal railing railing.
- Ø125 mm, diameter 360° swivel caster.
- Each wheel is individually braked.
- MDF lying surface.
- Disassembled head rails.
- Antibacterial liquid-proof hospital bed



DIMENSIONS



TECHNICAL SPECIFICATIONS

