

MYS-5310B

3 MOTOR ELECTRIC HOSPITAL BED



Meyosis[®]

www.meyosis.com

MYS-5310B



GENERAL SPECIFICATIONS

- Backrest, legrest and knee movements, height adjustment, vascular positions are controlled by 3 motors
- Positions are controlled by hand control (on head side)
- Removable Polypropylene (PP) plastic head and foot ends with colour panels
- Foldable aluminum side rails
- 4 independent and detachable ABS lying sections
- Antibacterial electrostatic powder coated steel frame, open chassis
- ABS plastic covered chasis
- Ø125 mm, 360° swivel wheels with cross brake system
- Accessory adaptors on the corners, 1 height adjustable IV Pole
- Urine bag holders and belt hooks on both sides
- Protective plastic bumpers on corners
- 26 Density foam mattress with antibacterial, waterproof and zippered cover

OPERATION ANGLES



Vascular
0° - 30°



Knee Movement



Yükseklik




Back Movement
0° - 70°



Leg Movement
0° - 30°

DIMENSIONS



TECHNICAL SPECIFICATIONS



100 kg



Max. 200 kg



3 Motors



110-220V
50-60hz