

MYS-540N

PEDIATRIC BED WITH 4 MOTORS




The logo for Meyosis, featuring a stylized blue waveform icon followed by the word "Meyosis" in a bold, blue, sans-serif font.

www.meyosis.com

MYS-540N

GENERAL SPECIFICATIONS

- Electrically adjustable height, backrest, knee break, Trendelenburg, and reverse Trendelenburg positions
- Electric Fowler position function
- Wired control panel
- 4-section sleeping platform
- 4-piece side rails made of easy-to-clean ABS material
- Removable ABS head and foot boards
- Casters with individual brake
- Electrostatic powder-coated metal frame
- IV pole
- Shock-absorbing corner bumpers
- Easily removable head and foot boards with locking/unlocking system
- Drop-down side rails with slow descent system using gas springs or equivalent systems



OPERATION ANGLES



Vascular
0° - 30°



Knee Movement



Back Movement
0° - 70°



Leg Movement
0° - 30°



Trendelenburg
0° - 12°




Height



Cardiac

DIMENSIONS



Meyosis®

TECHNICAL SPECIFICATIONS



130 kg



Max. 180 kg



4 Motors



110-220V
50-60hz