

# MYS-5410NE

4 MOTORS ICU HOSPITAL BED



Meyosis®

[www.meyosis.com](http://www.meyosis.com)

# MYS-5410NE



- ABS plastic covered chassis
- Ø125 mm, 360° swivel wheels with central break system and directional lock
- Accessory adaptors on the corners, 1 height adjustable IV Pole
- Urine bag holders and belt hooks on both sides
- Protective plastic bumpers on corners
- Alternative color options for side rails
- 26 Density foam mattress with antibacterial, waterproof and zippered cover

## OPTIONAL FEATURES

- Battery
- Manual CPR
- X-RAY
- Scale
- Side Control Panel

## GENERAL SPECIFICATIONS

- Backrest, legrest and knee movements, height adjustment, vascular, trendelenburg and reverse trendelenburg positions are controlled by 4 DC motors
- Positions are controlled by hand control (on head side)
- Electric CPR and function locking features on attendant control panel
- The dual auto-regression system, reduces the pressure on the abdomen and the friction on the pelvis area
- Removable Polypropylene (PP) plastic head and foot ends
- Split PP plastic side rails with soft drop realise mechanism and angle indicator
- 4 independent and detachable ABS lying sections
- Antibacterial electrostatic powder coated steel frame



### Sliding Backrest

The sliding backrest protects the patient and prevents uncomfortable forces during position adjustments.

### Position Safety

Side rails provide support for patients. The head-end grip bars make bed maneuvering easier.

### Collision Prevention

Wall bumper cones protect the bed and room furniture from damage.

### Control Panel

Positions can be adjusted electronically via the buttons on the control panel.

### Hygienic Mattress Bases

4-piece easy-to-clean ABS plastic mattress surface.

### Height Adjustment Range

The 40 cm – 80 cm height range helps prevent patient falls and ensures a comfortable working environment for healthcare staff.

### Central Brake

360° rotatable central brake system with directional lock castors.



## OPERATION ANGLES



Vascular  
0° - 30°



Knee Movement



Back  
Movement  
0° - 70°



Leg Movement  
0° - 30°



Trendelenburg  
0° - 17°




Height



Cardiac

## DIMENSIONS



## TECHNICAL SPECIFICATONS



125 kg



Max. 250 kg



4 Motors



110-220V  
50-60hz

