

MYS-650

EMERGENCY STRECHER




The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font. A registered trademark symbol (®) is located at the top right of the word.

www.meyosis.com

MYS-650



DIMENSIONS



TECHNICAL SPECIFICATIONS



GENERAL SPECIFICATIONS

- Backrest adjustable with gas spring, leg section adjustable with gear mechanism.
- Height, Trendelenburg and reverse Trendelenburg positions adjustable via hydraulic pedal pumps located on both sides.
- Full-length polyurethane side rails.
- Foldable push handles at the head and foot ends (metal).
- Full-length X-RAY cassette holder.
- Oxygen cylinder holder.
- 1 IV pole.
- Ø200 mm, 360° swivel castors with central locking and directional lock.
- Plastic corner bumpers on all 4 corners.
- Compact lying surface.
- Undercarriage covered with ABS plastic, metal parts with electrostatic powder coating.
- Waterproof mattress.

OPERATION ANGLES

