

# MYS-710T

OVER BED TABLE WITH ABS TOP  
AND TILT MECHANISM



The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font. A registered trademark symbol (®) is located at the top right of the word.

[www.meyosis.com](http://www.meyosis.com)

# MYS-710T



## GENERAL SPECIFICATIONS

- Whole ABS body
- ABS top with tilt mechanism
- Height adjustable with hydraulic
- 360 degrees swivel casters.

## DIMENSIONS

- |                  |           |
|------------------|-----------|
| — Width          | 42 cm     |
| — Length         | 71 cm     |
| — Min-Max Height | 82-102 cm |

