

MYS-718

MECHANICAL BLOOD DONOR CHAIR



The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font. A registered trademark symbol (®) is located at the top right of the word.

www.meyosis.com

MYS-718



GENERAL SPECIFICATIONS

- The backrest is adjusted with a shock absorber.
- The leg section can move together with the backrest or independently.
- A telescopic IV pole is included.
- The armrests move up and down by rotating.
- All metal parts are coated with electrostatic powder paint.
- There are 2 brake wheels available.

MAXIMUM LOAD CAPACITY



DIMENSIONS

Chair Position

- | | |
|----------|--------|
| — Width | 50 cm |
| — Depth | 80 cm |
| — Height | 120 cm |

Stretcher Position

- | | |
|----------|-------|
| — Width | 50 cm |
| — Length | 180cm |
| — Height | 55 cm |

