

MYS-728A

MEDICATION AND ANESTHESIA CART



Meyosis®

www.meyosis.com

MYS-728A



TECHNICAL SPECIFICATIONS

- The entire surface of the cart is made from hygienic material that is easy to clean.
- All surfaces are made of PP (polypropylene) material that is resistant to breakage and bending.
- Steel bushings are used at the wheel connection points.
- There is a foldable medication preparation tray on the side of the cart.
- On the right side of the cart, there is a removable and retractable writing board.

ERGONOMIC DESIGN

- The cart features:
 - 2 small drawers
 - 3 medium drawers for storing essential medical items.
- Drawer rails are spring-loaded and close automatically with a gentle push.
- Each drawer has a centrally located ergonomic handle integrated into the drawer front.
- The locking mechanism is centralized and can lock all drawers at once

FUNCTIONAL DRAWERS

- 1st and 2nd drawers: 16 compartments each
- 3rd, 4th, and 5th drawers: 9 compartments each
- Drawer interiors can be customized by the user without the need for any additional parts, according to preference and ease of use.
- Cart dimensions (depth x width): 58.75 cm x 75.4 cm
- Height from the ground: 86.25 cm

ACCESSORIES

- The cart includes 1 waste bin and 1 sharps container.
- There are 10 forward-tilting medication bins