

MYS-728N

EMERGENCY CART



The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font.

www.meyosis.com

MYS-728N



TECHNICAL SPECIFICATIONS

- The entire surface of the cart is made from hygienic material that is easy to clean.
- All surfaces are made of PP (polypropylene) material resistant to breaking and bending.
- Steel bushings are used at the wheel connection points.
- There is a foldable medication preparation tray on the side of the cart.
- The cart includes 1 waste bin and 1 sharps container.
- Equipped with 1 oxygen cylinder holder.
- Includes 1 resuscitation board.
- Includes 1 defibrillator shelf.
- Includes 1 IV pole.



ERGONOMIC DESIGN & FUNCTIONAL DRAWERS

- Drawer interiors can be customized by the user without any additional parts, according to preference and ease of use.
- Cart dimensions (depth x width): 58.75 cm x 75.4 cm
- Height from the ground: 86.25 cm

ACCESSORIES

- The cart contains:
 - 3 small drawers
 - 1 medium drawer
 - 1 large drawer for storing essential medical supplies.
- Drawer rails are spring-loaded and close automatically with a gentle push.
- Each drawer has a central ergonomic handle integrated into the drawer front for easy opening and closing.