

# MYS-728P

DRESSING TROLLEY



Meyosis<sup>®</sup>

[www.meyosis.com](http://www.meyosis.com)



# MYS-728P



## TECHNICAL SPECIFICATIONS

- The entire surface of the trolley is made from hygienic material that is easy to clean.
- All surfaces are made of PP (polypropylene) material that is resistant to breaking and bending.
- Steel bushings are used at the wheel connection points

## ERGONOMIC DESIGN

- The trolley includes:
  - 2 small drawers
  - 3 medium drawers
- 1 large drawer for storing medical supplies.
- Drawer rails are spring-loaded and close automatically with a gentle push.
- Each drawer has an ergonomic handle integrated into the drawer front

## FUNCTIONAL DRAWERS

- 1st and 2nd drawers: 16 compartments each
- 3rd, 4th, and 5th drawers: 9 compartments each
- The drawer locking mechanism is centralized and allows locking all drawers at once.
- Drawer interiors can be customized by the user without the need for additional equipment, according to preference and ease of use.
- Dimensions (depth x width): 58.75 cm x 75.4 cm
- Height from ground: 86.25 cm



## ACCESSORIES

- The cart includes 1 trash bin and 1 waste bin.
- On the right side of the cart, there is a removable and retractable writing board.