

# MYS-931

HAEMODIALYSIS AND CHEMOTHERAPY CHAIR



The logo for Meyosis, featuring a stylized blue waveform icon followed by the word "Meyosis" in a bold, blue, sans-serif font.

[www.meyosis.com](http://www.meyosis.com)

# MYS-931



## DIMENSIONS

Length	210 cm
Back Length	67,5 cm
Min-Max Height	56-76 cm
Seat Height	52 cm
Foot Section Length	62,5 cm
Footrest Extension	24,5 cm

## TECHNICAL SPECIFICATIONS



174 kg




Max. 200 kg

## GENERAL SPECIFICATIONS

- 4 Motors • Smart Control • Ergonomic Comfort
- Backrest, leg, height, and Trendelenburg positions can be easily adjusted via hand remote control
- Emergency mode: One-touch head-down position and return
- Synchronized armrests: Move with the backrest, adjustable up-down and side-to-side
- Cramp-preventing comfort with manually adjustable leg support
- Antibacterial, waterproof, non-sweating artificial leather upholstery
- All outer surfaces are covered with easy-to-clean ABS plastic
- The frame is corrosion-resistant and coated with electrostatic powder paint baked at 200°C
- Compatible with 220/380V 50Hz mains, protected against voltage fluctuations
- Optional battery backup for uninterrupted use during power outages
- Central brake, IV pole, and head pillow are standard accessories

## OPERATION ANGLES



## LEATHER COLOUR OPTIONS

Navy Blue	Beige
White	Pastel Blue
Black	Turquoise
Burgundy	Grey