

MYS-950

HYDRAULIC STRECHER



The logo for Meyosis, featuring a stylized blue ECG line that forms the letter 'M' followed by the word 'eyosis' in a blue sans-serif font, with a registered trademark symbol (®) to the upper right.

www.meyosis.com


MYS-950



GENERAL SPECIFICATIONS

- Back side adjustment by gas spring
- Height adjustment by hydraulic pedals at both sides
- Trendelenburg position is controlled by the pedal at foot side
- Full length aluminium side rails
- Antibacterial electrostatic powder coated steel frame
- ABS plastic covered chasis
- Accessory adaptors on the corners, 1 height adjustable IV Pole
- Ø150 mm, 360° swivel wheels with central brake system and directional lock
- Waterproof mattress with antibacterial, waterproof and zippered cover
- Wired basket on lower chasis for patient's belonging.

DIMENSIONS



TECHNICAL SPECIFICATIONS



OPERATION ANGLES

